Cornell Center for Materials Research

Educational Programs

Who's Got The Flu?

Author:Nuenyan LamDate Created:August 13, 2010Subject:Life ScienceLevel:8 - 12Standards:CCSS: Reading Standards for Literacy in Science 3, 6-9Schedule:Two 50 minute periods

Objectives:

The human body uses active immunity to defend against infectious diseases as a second-line of defense mechanism.
 Enzyme-linked immunoabsorbant assay (ELISA) may be used as a diagnostic tool to identify Infectious diseases.
 Students will:

- Understand how the human body defends against infectious diseases.
- Understand how ELISA is used a diagnostic tool to diagnose infectious diseases.
- Use an ELISA simulation kit to diagnose the Avian Flu Virus in hypothetical patients.
- Understand how infectious diseases are transmitted.

Vocabulary:

ELISA Antibody Antigen Immunity Active Immunity Passive Immunity Enzyme Substrate

Materials:

For Each Group:

- ELISA Simulation Kit from Carolina Biological:
 - One microtiter plate
 - Eleven disposable pipettes
 - One copy of patient profiles of each disease being tested
 - Student Guide worksheet for each student
 - 300 uL diluted solutions of positive and negative controls
 - 1500 uL diluted solutions of patient samples A-F
 - Diluted solutions of antigen, secondary antibody, and chromogen (ALL are simulated samples)
 - Internet access

Safety:

- All solutions used in this activity are simulated, thus there is no risk of infection of any disease.
 Do not ingest any of the solutions
 - Do not ingest any of the solutions used during this activity.
 - Actively supervise students throughout the entire laboratory investigation.

Science Content for the Teacher:

The human body uses two methods to defend against infectious diseases. Our skin serves as our first line of defense. It serves as a barrier between the internal and external environments. Cilia present in our respiratory tract also prevent pathogens from invading our body. Coughing and sneezing are ways our bodies respond to foreign particles that might have entered our bodies. However, when bacteria or viruses invade our bodies, the second line of defense is stimulated to help us fight off these, possibly, deadly pathogens.

Our second line of defense, also called active immunity, is a process in which are our body produces antibodies that identifies a certain pathogen. Once identified, a cascade of reactions occurs to stimulate the production of white blood cells to fight off and destroy the foreign pathogen.

Viruses and bacteria are coated with specific proteins that only certain antibodies can recognize. The enzyme-linked immunoabsorbant assay (ELISA) is a diagnostic tool to detect the specific proteins that may be present in patient blood samples. Antigen proteins are first placed in the wells of a microtiter plate. Patient samples with antibodies are added. The antibody will contain sites that are specific to the antigen proteins will adhere to each other and form a complex. A secondary antibody (anti-antibody) is added. A special enzyme is attached to this secondary antibody. This antibody will attach to the patient sample. In the final step, a substrate is added to attach to the enzyme to produce a color change in the sample to indicate a presence of the virus or bacteria in the sample.

Preparation:

- Dilute original patient samples and antibody solutions to allow enough samples for each group.
- Provide students with basic background in how active immunity occur in the body.
- Equip each group with necessary copies of Student Guides, Patient Profiles, and Periodicals.
- Ensure each group has one microtiter plate and enough pipettes.
- Review Lab Safety Rules with students





Educational Programs

Classroom Procedure:

Engage (Time: 40 min)

Engage students with primary sources from the World Health Organization. Since the outbreak of the avian influenza virus in 1997, the WHO has been tracking the occurrence of each case around the world.

http://gamapserver.who.int/mapLibrary/Files/Maps/Global_H5N1inHumanCUMULATI VE_FIMS_20100831.png

Provide students with periodical from the New York Times related to the concerns of the Avian Flu Virus, <u>http://www.nytimes.com/2004/11/30/international/asia/30flu.html?_r=1</u>

Allow students to imagine that they are immunologists who just received patient samples from a local hospital. Doctors suspect these patients are infected with the Avian Flu Virus. Students must use the ELISA test determine which patients are infected with the virus.

Model the procedure to perform the ELISA test. Be sure to remind students to change pipettes every time samples are transferred to avoid contamination.

Explore (Time: 40 min)

Following the directions presented in the Student Guide, perform ELISA tests on patient A through F.

Draw a diagram of the microtiter plate before placing solutions in each well. Have students plan which solutions should be added in each well of the plate to avoid confusion.

Have students practice using the pipettes with plain water first. Perform positive and negative control tests.

Perform tests on samples A through F with three trials for each sample. Record observations on Student Guide. A purple color change indicates a positive test and a light green color change indicates a negative test for the virus.

Explain (Time: 30 min)

Compare the results of each group. Ask students which patients are infected with the Avian Flu Virus. Explain to students the function of each solution added to each well during the ELISA test. Link the antibody-antigen complex to the enzyme-substrate complex. The final step is used as an indicator of the presence of the virus. Show students the following animation to differentiate learning: http://www.sumanasinc.com/webcontent/animations/content/ELISA.html







Expand (Time: 10 min)

During this activity, students performed an ELISA test to identify the Avian Flu Virus. Ask students conduct research to answer the following questions:

- What other diseases can ELISA detect?
- Can ELISA detect non-communicable diseases, such as diabetes or cancer?

Assessment:

The following rubric can be used to assess students during each part of the activity. The term "expectations" here refers to the content, process and attitudinal goals for this activity. Evidence for understanding may be in the form of oral as well as written communication, both with the teacher as well as observed communication with other students. Specifics are listed in the table below.

- 1= exceeds expectations
- 2= meets expectations consistently
- 3= meets expectations occasionally
- 4= not meeting expectations

	Engage	Explore	Explain	Expand/Synthesis
1	Shows leadership in the	Completes work	Provides an in-depth	Shows understanding of
	discussion and link	accurately while	explanation of findings,	how ELISA test can be
	real-world data to	providing an	making good use of	used to detect various
	current events relating	explanation for	vocabulary terms.	diseases. Writes detailed
	to the avian flu virus	what is observed.	Writes a detailed	description of the process
		Works very well	description in science	in science journal.
		with group.	journal.	
2	Participates in the	Completes work	Provides clear	Shows understanding and
	discussion and shows	accurately and	explanation of findings.	can express a good
	an understanding of	works	Writes a good	description of the process
	real-world examples of	cooperatively	description of	of performing an ELISA
	infectious diseases	with group.	observations and	test.
	transmission.		procedure.	
3	Contributes to the	Works	Provides a limited	Shows some understanding
	brainstorm, but shows	cooperatively	explanation of findings	of the process of an ELISA
	little understanding of	with group, but	in science journal.	test. Records the process in
	real-world data linked	makes some	-	science journal.
	to infectious diseases	mistakes with the		
	transmission.	procedure.		
4	Does not participate in	Has trouble	Is not clear in	Does not understand how
-	brainstorm. Shows no	working with	explanation of findings.	ELISA tests are performed.
	understanding of how	group. Does little	Incomplete explanation	Incomplete description of
	real data linked to	to complete the	in science journal.	the process.
	infectious diseases	procedure.	5	*
	transmission	1		





Extension Activities:

- In this activity, students used an "indirect ELISA test" to detect an infectious disease. How is a "direct ELISA test" different? Have students research the difference between a direct and indirect ELISA test.
- Create a flip-book showing a positive or negative ELISA test.

Supplemental Information:

- More information on the Avian Flu Virus at www.who.int
- Periodicals on the Avian Flu Virus available at www.nytimes.com
- Animations on ELISA available at <u>http://www.sumanasinc.com/webcontent/animations/content/ELISA.html</u> and <u>http://www.biology.arizona.edu/IMMUNOLOGY/activities/elisa/technique.html</u>

Safety:

- All solutions used in this activity are simulated, thus there is no risk of infection of any disease.
- Review lab safety rules before this laboratory investigation.
- Do not ingest any of the solutions used during this activity.
- Actively supervise students throughout the entire laboratory investigation.

Acknowledgments:

Dr. Antje Baeumner, Dr. Katie Edwards, Kevin Dilley, Nev Singhota, NSF



